

MiSight[®] 1 day

for daily wear

Product Reference Guide

Brand	MiSight [®] 1 day
Material/ H ₂ O Content	omafilcon A / 60%
Base Curve (mm)	8.7
Dia (mm)	14.2
Sphere Power	-0.50D to -7.00D (0.50D steps after -6.00D) No plano
FDA Group	2
Wearing Schedule	Daily Disposable
Dk/t (at -3.00D)	25
Revenue Carton Size	90-pack
Trials Labeled As/Pack Size	Myopia Management ActivControl™ Technology / 10pk

Features

- Visibility tint
- Designed for myopia control



CooperVision[®]